

Product datasheet for **RC224538L4V**

TXNDC11 (NM_015914) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	TXNDC11
Synonyms:	EFPI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015914
ORF Size:	2874 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224538).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_015914.5 , NP_056998.4
RefSeq Size:	3047 bp
RefSeq ORF:	2877 bp
Locus ID:	51061
UniProt ID:	Q6PKC3
Cytogenetics:	16p13.13



Domains: thioered

Protein Families: Druggable Genome

MW: 107.3 kDa

Gene Summary: May act as a redox regulator involved in DUOX proteins folding. The interaction with DUOX1 and DUOX2 suggest that it belongs to a multiprotein complex constituting the thyroid H₂O₂ generating system. It is however not sufficient to assist DUOX1 and DUOX2 in H₂O₂ generation.[UniProtKB/Swiss-Prot Function]