

Product datasheet for **RC224542L3V**

B3GALNT1 (NM_003781) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | B3GALNT1 (NM_003781) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | B3GALNT1 |
| Synonyms: | B3GALT3; beta3Gal-T3; galT3; Gb4Cer; GLCT3; GLOB; P; P1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_003781 |
| ORF Size: | 993 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC224542). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003781.3 |
| RefSeq Size: | 3319 bp |
| RefSeq ORF: | 996 bp |
| Locus ID: | 8706 |
| UniProt ID: | O75752 |
| Cytogenetics: | 3q26.1 |
| Domains: | Galactosyl_T |
| Protein Families: | Transmembrane |



[View online »](#)

Protein Pathways: Glycosphingolipid biosynthesis - globo series, Metabolic pathways

MW: 39.5 kDa

Gene Summary: This gene is a member of the beta-1,3-galactosyltransferase (beta3GalT) gene family. This family encodes type II membrane-bound glycoproteins with diverse enzymatic functions using different donor substrates (UDP-galactose and UDP-N-acetylglucosamine) and different acceptor sugars (N-acetylglucosamine, galactose, N-acetylgalactosamine). The beta3GalT genes are distantly related to the Drosophila Brainiac gene and have the protein coding sequence contained in a single exon. The beta3GalT proteins also contain conserved sequences not found in the beta4GalT or alpha3GalT proteins. The carbohydrate chains synthesized by these enzymes are designated as type 1, whereas beta4GalT enzymes synthesize type 2 carbohydrate chains. The ratio of type 1:type 2 chains changes during embryogenesis. By sequence similarity, the beta3GalT genes fall into at least two groups: beta3GalT4 and 4 other beta3GalT genes (beta3GalT1-3, beta3GalT5). The encoded protein of this gene does not use N-acetylglucosamine as an acceptor sugar at all. [provided by RefSeq, Mar 2017]