

Product datasheet for **RC224553L3V**

ALKBH6 (NM_032878) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ALKBH6 (NM_032878) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ALKBH6
Synonyms:	ABH6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032878
ORF Size:	798 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224553).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032878.3 , NP_116267.3
RefSeq Size:	980 bp
RefSeq ORF:	801 bp
Locus ID:	84964
UniProt ID:	Q3KRA9
Cytogenetics:	19q13.12
Domains:	2OG-Fell_Oxy
MW:	29.1 kDa



[View online »](#)

Gene Summary: Probable dioxygenase that requires molecular oxygen, alpha-ketoglutarate and iron.
[UniProtKB/Swiss-Prot Function]