

Product datasheet for **RC224577**

PDE9A (NM_001001577) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDE9A (NM_001001577) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PDE9A
Synonyms:	HSPDE9A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC224577 representing NM_001001577
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGACGCATTAGAAGAACACGACCATCTCCCTGCTGACCACCGACGACGCCATGGTCTCCATCGACC
 CCACCATGCCCGGAATTCAGAACGAGCATTCAAATCAATGAACTGAAAGCTGAAGTTGCAAACTACTT
 GGCTGTCTAGAGAAACGCGTGGAAATTGGAAGGACTAAAAGTGGTGGAGATTGAGAAATGCAAGAGTGAC
 ATTAAGAAGATGAGGGAGGAGCTGGCGGCCAGAAGCAGCAGGACCAACTGCCCTGTAAGTACAGTTTTT
 TGGATAACCACAAGAAGTTGACTCCTCGACGCGATGTTCCCACTTACCCCAAGTACCTGCTCTCCAGA
 GACCATCGAGGCCCTGCGGAAGCCGACCTTTGACGTCTGGCTTTGGGAGCCCAATGAGATGCTGAGCTGC
 CTGGAGCACATGTACCACGACCTCGGGTGGTCAAGGACTTCCAGCATCAACCCTGTACCCTCAGGAGGT
 GGCTGTTCTGCGTCCACGACAACACTACAGAAACAACCCCTTCCACAACCTCCGGCACTGCTTCTGCGTGGC
 CCAGATGATGTACAGCATGGTCTGGCTCTGCAGTCTCCAGGAGAAGTTCTCACAACGGATATCCTGATC
 CTAAATGACAGCGCCATCTGCCACGATCTGGACCATCCCGCTACAACAACACGTACCAGATCAATGCC
 GCACAGAGCTGGCGGTCCGCTACAATGACATCTCACCGCTGGAGAACCACCACTGCGCCGTGGCCTTCCA
 GATCCTCGCCGAGCCTGAGTGAACATCTTCTCCAACATCCCACCTGATGGGTTCAAGCAGATCCGACAG
 GGAATGATCACATTAATCTTGGCCACTGACATGGCAAGACATGCAGAAATATGGATTCTTTCAAAGAGA
 AAATGGAGAATTTGACTACAGCAACGAGGAGCACATGACCCTGCTGAAGATGATTTTGATAAAATGCTG
 TGATATCTTAACGAGGTCGGTCCAATGGAAGTCGAGAGCCTTGGTGGACTGTTTATTAGAGGAATAT
 TTTATGCAGAGCGACCGTGAAGAAGTCAGAAGCCCTTCTGTGGCACCGTTCATGGACCGAGACAAAGTGA
 CCAAGGCCACAGCCAGATTGGGTTTCATCAAGTTTGTCTGATCCCAATGTTTGAACAGTGACCAAGTGA
 CTTCCCATGGTTGAGGAGATCATGCTGCAGCCACTTTGGGAATCCCGAGATCGCTACGAGGAGCTGAAG
 CGGATAGATGACCCATGAAAGAGTTACAGAAGAAGACTGACAGCTTGACGTCTGGGGCCACCGAGAAGT
 CCAGAGAGAGAAGCAGAGATGTAAAAACAGTGAAGGAGACTGTGCC

AC**CGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA**
 TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC224577 representing NM_001001577
 Red=Cloning site Green=Tags(s)

MDA~~FR~~NTTISLLTTDDAMVSI~~DP~~TMPANSE~~RAF~~KINELKAEVANHLAVLEK~~RVE~~LEGLKVVEIEKCKSD
 IKKMREELAARSSRTNCPCKYSFLDNHKKLTPRRDVPTYPKYLLSPETIEALRKP~~TF~~DVWLWEPNEMLSC
 LEHMYHDLGLVRDFSINPVTLRRWLFVHDNYRNNPFHNRHCFV~~QA~~MMYSMVWLCSLQEKFSQTDILI
 LMTAAICHDLDPGYNN~~TY~~QINARTELAVRYNDISPLENHHCAVAFQILAEPECNIFSNIP~~PD~~GFKQIRQ
 GMITLILATDMARHAEIMDSFKEKMF~~ND~~YSNEEHMTLLKMILIKCCDISNEVRPMEVAEPWVDCLLEEY
 FMQSDREKSEGLPVAPFMDRDKVTKATAQIGFIK~~FV~~LIPMFETVTKLFP~~VEE~~IMLQPLWESRDRYEELK
 RIDDAMKELQKKTDSLTS~~GA~~TEKSRERSR~~VD~~KNSEGDCA

TRRLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk8061_d09.zip

Restriction Sites:

Sgfl-NotI

Cloning Scheme:



ACCN: NM_001001577

ORF Size: 1377 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001577.1](#), [NP_001001577.1](#)

RefSeq Size: 1753 bp

RefSeq ORF: 1380 bp

Locus ID: 5152

UniProt ID: [O76083](#)

Cytogenetics: 21q22.3


Protein Families: Druggable Genome

Protein Pathways: Progesterone-mediated oocyte maturation, Purine metabolism


MW: 53.3 kDa

Gene Summary: The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC224577



Western blot validation of overexpression lysate (Cat# [LY424394]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC224577 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).