

## Product datasheet for **RC224582L1V**

### Carboxypeptidase M (CPM) (NM\_001874) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Carboxypeptidase M (CPM) (NM_001874) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Carboxypeptidase M
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001874
ORF Size:	1329 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224582).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001874.3</a>
RefSeq Size:	6683 bp
RefSeq ORF:	1332 bp
Locus ID:	1368
UniProt ID:	<a href="#">P14384</a>
Cytogenetics:	12q15
Domains:	Zn_carbOpept
Protein Families:	Druggable Genome, Protease
MW:	50.51 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a membrane-bound arginine/lysine carboxypeptidase. Its expression is associated with monocyte to macrophage differentiation. This encoded protein contains hydrophobic regions at the amino and carboxy termini and has 6 potential asparagine-linked glycosylation sites. The active site residues of carboxypeptidases A and B are conserved in this protein. Three alternatively spliced transcript variants encoding the same protein have been described for this gene. [provided by RefSeq, Jul 2008]