

Product datasheet for **RC224583**

PPP2R2B (NM_181676) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP2R2B (NM_181676) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PPP2R2B
Synonyms:	B55BETA; PP2AB55BETA; PP2ABBETA; PP2APR55B; PP2APR55BETA; PR2AB55BETA; PR2ABBETA; PR2APR55BETA; PR52B; PR55-BETA; PR55BETA; SCA12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC224583 representing NM_181676
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAAATGCTTCTCTCGTTACCTGCCTTACATCTTCAGACCTCCGAACACCATCCTTTCTCCAGCTGCC
 ACACAGAAGCTGACATTATCTCTACGGTAGAATTC AACCCACACGGGAGAATTA CTAGCGACAGGGGACAA
 GGGGGTTCGGGTTGTAATATTTCAACGAGAGCAGGAGAGTAAAAATCAGGTTTCATCGTAGGGGTGAATAC
 AATGTTTACAGCACATTCAGAGCCATGAACCCGAGTTCGATTACCTGAAGAGTTTAGAAAATAGAAGAAA
 AAATCAATAAAAATAAGATGGCTCCCCAGCAGAAATGCAGCTTACTTTCTTCTGTCTACTAATGATAAAAC
 TGTGAAGCTGTGAAAAGTCAGCGAGCGTGATAAGAGGCCAGAAGGCTACAATCTGAAAGATGAGGAGGGC
 CGGCTCCGGGATCCTGCCACCATCACACCCTGCGGGTGCCTGTCCTGAGACCCATGGACCTGATGGTGG
 AGGCCACCCACGAAGAGTATTTGCCAACGCACACATATCACATCACTCCATATCTGTCAACAGCGA
 CTATGAAACCTACATGTCGCTGATGACCTGAGGATTAACCTATGGAACCTTGAATAACCAATCAAAGT
 TTTAATATTGTGGACATTAAGCCAGCCAACATGGAGGAGCTCACGGAGGTGATCACAGCAGCCGAGTTCC
 ACCCCCATCATTGCAACACCTTCGTGTACAGCAGCAGCAAAGGGACAATCCGGCTGTGTGACATGCGGGC
 ATCTGCCCTGTGTGACAGGCACACCAAATTTTTGAAGAGCCGGAAGATCCAAGCAACAGATCATTTTTTC
 TCTGAAATATCTCTTCGATTTCCGATGTGAAGTTCAGCCACAGTGGGAGGTATATCATGACCAGGGACT
 ACTTGACCGTCAAAGTCTGGGATCTCAACATGGAAAACCGCCCCATCGAGACTTACCAGGTTTCATGACTA
 CCTCCGCAGCAAGCTGTGTTCCCTCTATGAAAATGACTGCATTTTTGATAAATTTGAGTGTGTGGAAT
 GGGTCAGACAGTGTATCATGACAGGCTCTACAACAACCTTTCAGGATGTTTCGACAGAAACACCAAGC
 GTGATGTGACCCTTGAGGCTTCGAGGGAAAACAGCAAGCCCGGGCTATCCTCAAACCCGAAAAGTGTG
 TGTGGGGGCAAGCGGAGAAAAGACGAGATCAGTGTGACAGTCTGGACTTTAGCAAAAAGATCTTGCAT
 ACAGCTTGGCATCCTTCAGAAAATATTATAGCAGTGGCGCTACAATAACCTATATATATTCCAGGACA
 AGGTTAAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC224583 representing NM_181676
 Red=Cloning site Green=Tags(s)

MKCFSRYPYIFRPPNTILSSSCHTEADIISTVEFNHTGELLATGDKGGRVVIFQREQESKNQVHRRGEY
 NVYSTFQSHEPEFDYLSLEIEEKINKIRWLPQQNAAYFLLSTNDKTVKLVKVSERDKRPEGYNLKDEEG
 RLRDPATITTLRVPVLRPMDLMVEATPRRVFANAHTYHINSISVNSDYETYMSADDLRINLWNFEITNQS
 FNIVDIK PANMEELTEVITAAEFHPHCNTFVYSSSKGTIRLCDMRASALCDRHTKFFEEPDPNSRFF
 SEI ISSISDVKFSHSGRYIMTRDYLTVKVWDLNMENRPIETYQVHDYLRSLKCSLYENDCIFDKFECVWN
 GSDSVIMTGSYNNFRMFDNRNTRDVTLEASRENSKPRAILKPRKVCVGGKRRKDEISVDSLDFSKKILH
 TAWHPSENIIVAATNNLYIFQDKVN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8060_f01.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_181676

ORF Size: 1338 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181676.2](#)

RefSeq Size: 1978 bp

RefSeq ORF: 1341 bp

Locus ID: 5521

UniProt ID: [Q00005](#)

Cytogenetics: 5q32


Protein Families: Druggable Genome, Phosphatase

Protein Pathways: Tight junction


MW: 51.8 kDa

Gene Summary: The product of this gene belongs to the phosphatase 2 regulatory subunit B family. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes a beta isoform of the regulatory subunit B55 subfamily. Defects in this gene cause autosomal dominant spinocerebellar ataxia 12 (SCA12), a disease caused by degeneration of the cerebellum, sometimes involving the brainstem and spinal cord, and in resulting in poor coordination of speech and body movements. Multiple alternatively spliced variants, which encode different isoforms, have been identified for this gene. The 5' UTR of some of these variants includes a CAG trinucleotide repeat sequence (7-28 copies) that can be expanded to 55-78 copies in cases of SCA12. [provided by RefSeq, Jul 2016]

Product images:



Circular map for RC224583



Western blot validation of overexpression lysate (Cat# [LY405665]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC224583 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).