

Product datasheet for **RC224589L3V**

SCRT1 (NM_031309) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SCRT1
Synonyms:	SCRT; ZNF898
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031309
ORF Size:	1044 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224589).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_031309.3
RefSeq Size:	3552 bp
RefSeq ORF:	1047 bp
Locus ID:	83482
UniProt ID:	Q9BWW7
Cytogenetics:	8q24.3



MW: 35.4 kDa

Gene Summary: This gene encodes a C2H2-type zinc finger transcriptional repressor that binds to E-box motifs. The encoded protein may promote neural differentiation and may be involved in cancers with neuroendocrine features. [provided by RefSeq, Jul 2013]