

## Product datasheet for **RC224605L4V**

### SUN1 (NM\_025154) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SUN1 (NM_025154) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SUN1
Synonyms:	UNC84A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025154
ORF Size:	2106 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224605).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025154.5</a>
RefSeq Size:	4047 bp
RefSeq ORF:	2109 bp
Locus ID:	23353
UniProt ID:	<a href="#">O94901</a>
Cytogenetics:	7p22.3
Protein Families:	Transmembrane
MW:	77.7 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the unc-84 homolog family and encodes a nuclear envelope protein with an Unc84 (SUN) domain. The protein is involved in nuclear anchorage and migration. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2019]