

Product datasheet for **RC224611A1V**

Human PSMB8 (NM_148919) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PSMB8 (NM_148919) AAV Particle
Tag:	Myc-DDK
Symbol:	PSMB8
Synonyms:	ALDD; D6S216; D6S216E; JMP; LMP7; NKJO; PRAAS1; PSMB5i; RING10
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC224611 representing NM_148919 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGCTACTAGATGTATGCGGAGCCCCCGAGGGCAGCGGCCGGAATCGGCTCTCCCGTTGCGGGAA
GCGGGCGTCGCTCGGACCCAGGACACTACAGTTTCTCTATGCGATCTCCAGAGCTCGCTTTACCCCGGG
AATGCAGCCACAGAATTCTCCAGTCCCTGGGTGGGACGGAGAAAGGAACGTTTCAGATTGAGATGGCC
CATGGCACCACCGCTCGCCTTCAAGTTCAGCATGGAGTGATTGCAGCAGTGGATTCTCGGSCCTCAG
CTGGGTCTACATTAGTGCCTTACGGGTGAACAAGGTGATTGAGATTAACCTTACCTGCTTGGCACCAT
GTCTGGCTGTGCAGCAGACTGTCACTACTGGGAGCGCCTGCTGGCCAAGGAATGCAGGCTGTACTATCTG
CGAAATGGAGAACGTATTTCACTGTGCGCAGCCTCCAAGCTGCTGTCCAACATGATGTGCCAGTACCGGG
GCATGGGCCTCTCTATGGGCAGTATGATCTGTGGCTGGGATAAGAAGGGTCTGGACTCTACTACGTGGA
TGAACATGGGACTCGGCTCTCAGGAAATATGTTCTCCACGGGTAGTGGGAACACTTATGCCTACGGGGTC
ATGGACAGTGGCTATCGGCCTAATCTTAGCCCTGAAGAGGCCTATGACCTTGGCCGAGGGCTATTGCTT
ATGCCACTCACAGAGACAGCTATTCTGGAGGCGTTGTCAATATGTACCACATGAAGGAAGATGGTTGGT
GAAAGTAGAAAGTACAGATGTCACTGACCTGCTGCACCACTACCGGGAAGCCAATCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC224611 representing NM_148919 Red=Cloning site Green=Tags(s) MALLDVCGAPRGQRPESALPVAGSGRRSDPGHYSFSMRSPALPRGMQPTTEFFQSLGGDGERNVQIEMA HGTTLAFKFQHGVIAAVDSRASAGSYISALRVNKVIEINPYLLGTMSGCAADCQYWERLLAKECRLYYL RNGERISVSAASKLLSNMMCQYRGMGLSMGSMICGWDKKGPLYVDEHGTRLSGNMFSTGSGNTYAYGV MDSGYRPNLSPEEAYDLGRRAIAYATHRDSYSGGVNMYHMKEDGWVKVESTDVSDLLHQYREANQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_148919
ORF Size:	828 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_148919.3</u>
RefSeq Size:	1135 bp
RefSeq ORF:	831 bp
Locus ID:	5696
UniProt ID:	<u>P28062</u>
Cytogenetics:	6p21.32
MW:	30.2 kDa