

Product datasheet for **RC224637L3V**

CCDC62 (NM_032573) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	CCDC62
Synonyms:	aaa; CT109; ERAP75; FLJ25432; FLJ40344; TSP-NY
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032573
ORF Size:	2046 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224637).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_032573.3](#), [NP_115962.2](#)

RefSeq Size: 2481 bp

RefSeq ORF: 2048 bp

Locus ID: 84660

Cytogenetics: 12q24.31

MW: 77.3 kDa



Gene Summary:

Nuclear receptor coactivator that can enhance preferentially estrogen receptors ESR1 and ESR2 transactivation. Modulates also progesterone/PGR, glucocorticoid/NR3C1 and androgen/AR receptors transactivation, although at lower level; little effect on vitamin D receptor/VDR.[UniProtKB/Swiss-Prot Function]