

Product datasheet for **RC224673L4V**

ZNF295 (ZBTB21) (NM_001098403) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF295 (ZBTB21) (NM_001098403) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF295
Synonyms:	ZNF295
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001098403
ORF Size:	2595 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224673).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001098403.1
RefSeq Size:	6876 bp
RefSeq ORF:	2598 bp
Locus ID:	49854
UniProt ID:	Q9ULJ3
Cytogenetics:	21q22.3
Protein Families:	Transcription Factors
MW:	95.6 kDa


[View online »](#)

Gene Summary: Acts as a transcription repressor.[UniProtKB/Swiss-Prot Function]