

Product datasheet for RC224716L3

PLEKHG3 (NM_015549) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PLEKHG3 (NM_015549) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	PLEKHG3
Synonyms:	ARHGEF43; KIAA0599
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224716).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 26030
Cytogenetics: 14q23.3
MW: 128.1 kDa