

Product datasheet for **RC224744L1V**

TIMM50 (NM_001001563) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TIMM50 (NM_001001563) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TIMM50
Synonyms:	MGCA9; TIM50; TIM50L
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001001563
ORF Size:	1368 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224744).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001001563.1 , NP_001001563.1
RefSeq Size:	2572 bp
RefSeq ORF:	1062 bp
Locus ID:	92609
UniProt ID:	Q3ZCQ8
Cytogenetics:	19q13.2
MW:	50.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a subunit of the TIM23 inner mitochondrial membrane translocase complex. The encoded protein functions as the receptor subunit that recognizes the mitochondrial targeting signal, or presequence, on protein cargo that is destined for the mitochondrial inner membrane and matrix. This protein may also play a role in maintaining the membrane permeability barrier, and knockdown of this gene in human cells results in the release of cytochrome c and apoptosis. [provided by RefSeq, Jul 2016]