

Product datasheet for RC224784A1V

Human IL36 gamma (IL36G) (NM_019618) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human IL36 gamma (IL36G) (NM_019618) AAV Particle
Tag: Myc-DDK
Symbol: IL36 gamma
Synonyms: IL-1F9; IL-1H1; IL-1RP2; IL1E; IL1F9; IL1H1; IL1RP2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC224784 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGAGAGGCACTCCAGGAGACGCTGATGGTGGAGGAAGGGCCGTCTATCAATCAATGTGTAAACCTATTA
 CTGGGACTATTAATGATTTGAATCAGCAAGTGTGGACCCTTCAGGGTCAGAACCTTGTCAGTTCCACG
 AAGTGACAGTGTGACCCAGTCACTGTTGCTGTTATCACATGCAAGTATCCAGAGGCTCTTGAGCAAGGC
 AGAGGGGATCCCATTTATTTGGGAATCCAGAATCCAGAAATGTGTTTGTATTGTGAGAAGGTTGGAGAAC
 AGCCCATTTGCAGCTAAAAGAGCAGAAGATCATGGATCTGTATGGCCAACCCGAGCCCGTGAAACCCCTT
 CCTTTTCTACCGTGCCAAGACTGGTAGGACCTCCACCCTTGAGTCTGTGGCCTTCCCGGACTGGTTCATT
 GCCTCCTCCAAGAGAGACCAGCCCATCATTCTGACTTCAGAACTTGGGAAGTCATACAACACTGCCTTTG
 AATTAATATAAATGAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224784 protein sequence
 Red=Cloning site Green=Tags(s)

MRGTPGDADGGGRAVYQSMCKPITGTINDLNQQVWTLQGQNLVAVPRSDSVPTVAVITCKYPEALEQG
 RGDPIYLGIQNPENCLYCEKVGEQPTLQLKEQKIMDLYGQPEPVKPFIFYRAKTGRTSTLESVAFPDWFI
 ASSKRDPPIILTSELGKSYNTAFELNIND

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_019618
ORF Size:	507 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_019618.2</u>
RefSeq Size:	1212 bp
RefSeq ORF:	510 bp
Locus ID:	56300
UniProt ID:	<u>Q9NZH8</u>
Cytogenetics:	2q14.1
MW:	18.7 kDa