

Product datasheet for **RC224798L4V**

CEP63 (NM_001042384) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	CEP63
Synonyms:	SCKL6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001042384
ORF Size:	1425 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224798).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001042384.1
RefSeq Size:	1906 bp
RefSeq ORF:	1428 bp
Locus ID:	80254
UniProt ID:	Q96MT8
Cytogenetics:	3q22.2



MW: 55.4 kDa

Gene Summary: This gene encodes a protein with six coiled-coil domains. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells. Several alternatively spliced transcript variants have been found, but their biological validity has not been determined. [provided by RefSeq, Jul 2008]