

## Product datasheet for **RC224807L3V**

### HN1 (JPT1) (NM\_001002032) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | HN1 (JPT1) (NM_001002032) Human Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                               |
| Symbol:                   | JPT1                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Synonyms:                 | ARM2; HN1; HN1A                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                                                                                                                                                                                                                                                                                                                                                                                                               |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_001002032                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| ORF Size:                 | 543 bp                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC224807).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_001002032.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                     |
| RefSeq Size:              | 1582 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 546 bp                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Locus ID:                 | 51155                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| UniProt ID:               | <a href="#">Q9UK76</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:             | 17q25.1                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| MW:                       | 19.7 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                           |


[View online »](#)

**Gene Summary:**

Modulates negatively AKT-mediated GSK3B signaling (PubMed:21323578, PubMed:22155408). Induces CTNNB1 'Ser-33' phosphorylation and degradation through the suppression of the inhibitory 'Ser-9' phosphorylation of GSK3B, which represses the function of the APC:CTNNB1:GSK3B complex and the interaction with CDH1/E-cadherin in adherent junctions (PubMed:25169422). Plays a role in the regulation of cell cycle and cell adhesion (PubMed:25169422, PubMed:25450365). Has an inhibitory role on AR-signaling pathway through the induction of receptor proteosomal degradation (PubMed:22155408). [UniProtKB/Swiss-Prot Function]