

Product datasheet for **RC224808L1V**

SCPI (SYCP1) (NM_003176) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SCPI
Synonyms:	CT8; HOM-TES-14; SCP-1; SCPI
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003176
ORF Size:	2928 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224808).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_003176.2
RefSeq Size:	3529 bp
RefSeq ORF:	2931 bp
Locus ID:	6847
UniProt ID:	Q15431
Cytogenetics:	1p13.2



MW: 114.2 kDa

Gene Summary: Major component of the transverse filaments of synaptonemal complexes, formed between homologous chromosomes during meiotic prophase. Required for normal assembly of the central element of the synaptonemal complexes. Required for normal centromere pairing during meiosis. Required for normal meiotic chromosome synapsis during oocyte and spermatocyte development and for normal male and female fertility.[UniProtKB/Swiss-Prot Function]