

Product datasheet for **RC224829L1V**

SNX30 (NM_001012994) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SNX30
Synonyms:	ATG24A
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001012994
ORF Size:	1311 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224829).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001012994.1
RefSeq Size:	7616 bp
RefSeq ORF:	1314 bp
Locus ID:	401548
UniProt ID:	Q5VWJ9
Cytogenetics:	9q32



MW: 49.5 kDa

Gene Summary: May be involved in several stages of intracellular trafficking.[UniProtKB/Swiss-Prot Function]