

Product datasheet for **RC224847L1V**

NRDE2 (NM_017970) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NRDE2 (NM_017970) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NRDE2
Synonyms:	C14orf102
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_017970
ORF Size:	3492 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224847).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017970.2
RefSeq Size:	3824 bp
RefSeq ORF:	3495 bp
Locus ID:	55051
UniProt ID:	Q9H7Z3
Cytogenetics:	14q32.11
MW:	132.5 kDa



[View online »](#)

Gene Summary:

Protein of the nuclear speckles that regulates RNA degradation and export from the nucleus through its interaction with MTREX an essential factor directing various RNAs to exosomal degradation (PubMed:30842217). Changes the conformation of MTREX, precluding its association with the nuclear exosome and interaction with proteins required for its function in RNA exosomal degradation (PubMed:30842217). Negatively regulates, for instance, the degradation of mRNAs and lncRNAs by inhibiting their MTREX-mediated recruitment to nuclear exosome (PubMed:30842217). By preventing the degradation of RNAs in the nucleus, it promotes their export to the cytoplasm (PubMed:30842217).[UniProtKB/Swiss-Prot Function]