

Product datasheet for **RC224861L3V**

CTBP2 (NM_001329) Human Tagged ORF Clone Lentiviral Particle


Product data:

Product Type:	Lentiviral Particles
Product Name:	CTBP2 (NM_001329) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTBP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001329
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224861).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001329.1
RefSeq Size:	2368 bp
RefSeq ORF:	1338 bp
Locus ID:	1488
UniProt ID:	P56545
Cytogenetics:	10q26.13
Domains:	2-Hacid_DH, 2-Hacid_DH_C
Protein Families:	Stem cell - Pluripotency, Stem cell relevant signaling - Wnt Signaling pathway



[View online »](#)

Protein Pathways:	Chronic myeloid leukemia, Notch signaling pathway, Pathways in cancer, Wnt signaling pathway
MW:	48.8 kDa
Gene Summary:	This gene produces alternative transcripts encoding two distinct proteins. One protein is a transcriptional repressor, while the other isoform is a major component of specialized synapses known as synaptic ribbons. Both proteins contain a NAD ⁺ binding domain similar to NAD ⁺ -dependent 2-hydroxyacid dehydrogenases. A portion of the 3' untranslated region was used to map this gene to chromosome 21q21.3; however, it was noted that similar loci elsewhere in the genome are likely. Blast analysis shows that this gene is present on chromosome 10. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Feb 2014]

Product images:

[RC224861L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC224861L3V particle to overexpress human CTBP2-Myc-DDK fusion protein.