

Product datasheet for **RC224863**

PDPN (NM_001006625) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PDPN (NM_001006625) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PDPN
Synonyms: AGGRUS; GP36; Gp38; GP40; HT1A-1; OTS8; PA2.26; T1A; T1A-2; T1A2; T11A
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC224863 representing NM_001006625
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCAGGTGCCGAAGATGATGTGGTGACTCCAGGAACCAGCGAAGACCGCTATAAGTCTGGCTTGACAA
CTCTGGTGCCAACAAGTGTCAACAGTGTAAACAGGCATTCGCATCGAGGATCTGCCAACTCAGAAAGCAC
AGTCCACGCGCAAGAACAAGTCCAAGCGCCACAGCCTCAAACGTGGCCACCAGTCACTCCACGGAGAAA
GTGGATGGAGACACAGACAACAGTTGAGAAAGATGGTTTGTCAACAGTGACCCGTGGTGAATCATAG
TTGGGGTCTTACTAGCCATCGGCTTCATTGGTGCAATCATCGTTGTGGTTATGCGAAAAATGTCGGGAAG
GCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224863 representing NM_001006625
Red=Cloning site Green=Tags(s)
MPGAEDDVVTPGTSEDRYKSLTTLVATSVNSVTGIRIEDLPTSESTVHAQEQSPSATASNVATSHSTEK
VDGDTQTTVEKDGLSTVTLVGIIVGVLLAIGFIGAIIVVVMRKMGRP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001006625

ORF Size: 354 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001006625.2
RefSeq Size:	2646 bp
RefSeq ORF:	357 bp
Locus ID:	10630
UniProt ID:	Q86YL7
Cytogenetics:	1p36.21
Protein Families:	Druggable Genome, Transmembrane
MW:	12.2 kDa
Gene Summary:	This gene encodes a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character. The homologous protein in other species has been described as a differentiation antigen and influenza-virus receptor. The specific function of this protein has not been determined but it has been proposed as a marker of lung injury. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]