

Product datasheet for **RC224876L4V**

ANKRD36B (NM_025190) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ANKRD36B
Synonyms:	KIAA1641
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025190
ORF Size:	4059 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224876).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_025190.3 , NP_079466.3
RefSeq Size:	5986 bp
RefSeq ORF:	4062 bp
Locus ID:	57730
UniProt ID:	Q8N2N9
Cytogenetics:	2q11.2



MW: 154 kDa