

Product datasheet for **RC224880L3V**

ST20 (NM_001100879) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ST20
Synonyms:	HCCS-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001100879
ORF Size:	237 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224880).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001100879.1](#), [NP_001094349.1](#)

RefSeq Size: 1025 bp

RefSeq ORF: 240 bp

Locus ID: 400410

Cytogenetics: 15q25.1

MW: 8.8 kDa



Gene Summary:

May act as a tumor suppressor. Promotes apoptosis of cancer cells.[UniProtKB/Swiss-Prot Function]