

Product datasheet for **RC224882A1V**

Human OR1D5 (NM_014566) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR1D5 (NM_014566) AAV Particle
Tag:	Myc-DDK
Symbol:	OR1D5
Synonyms:	C17orf2; OR1D4; OR17-2; OR17-30; OR17-31
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC224882 representing NM_014566 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATGGAGATAACCAGAGTGAGAACTCACAGTTCCTTCTCCTGGGGATCTCAGAGAGTCTGAGCAGC
AGCGGATCCTGTTTTGGATGTTCTGTCCATGTACCTGGTCACGGTGCTGGGAAATGTGCTCATCATCCT
GGCCATCAGCTCTGATTCCACCTGCACACCCCATGTACTTCTTCTGGCCAACTCTCCTTCACTGAC
CTCTTCTTTGTACCAACACAATCCCCAAGATGCTGGTGAACCTCCAGTCCAGAACAAGCCATCTCCT
ATGCAGGGTGTCTGACACAGCTCTACTTCTGGTCTCCTTGGTGACCTGGACAACCTCATCTGGCCGT
GATGGCGTATGATCGCTATGTGGCCACCTGCTGCCCCCTCCACTATGTCACAGCCATGAGCCCTGGGCTC
TGTGTCTTGTCTCTCCTTGTGTTGGGGCTGTCTGTTCTCTATGGCCTCCTCCTCACCTTCTCCTGA
CCAGGGTGACCTTCTGTGGGCTCGAGAGATCCACTACCTCTTCTGTGACATGTACATCCTGCTGTGGCT
GGCATGTTCCAACACCCACATCATTACACAGCGTTGATTGCCACTGGCTGCTTCATCTTCTCACCCCC
TTAGGGTTCATGACCACATCCTATGTACGTATTGTGAGAACCATCCTTCAAATGCCCTCGGCCTAAGA
AATAAAAACCTTCTCTACCTGTGCCTCCCATTTGGGTGTGGTCTCCCTCTTTATGGGACGCTTGCTAT
GGTGTACCTGCAGCCCTCCATACCTACTCCATGAAGGACTCAGTAGCCACAGTGATGTATGCTGTGCTG
ACACCTATGATGAACCTTTCTACAGGCTGAGGAACAAAGACATGCATGGGCTCCGGAAGAGTCC
TATGGAGACCCTTTCAGAGGCCTAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>RC224882 representing NM_014566 Red=Cloning site Green=Tags(s) MDGDNQSENSQFLLLGISESPEQQRILFWMFLSMYLVTVLGNVLIILAISSDSHLHTPMYFFLANLSFTD LFFVTNTIPKMLVNFQSQNKAI SYAGCLTQLYFLVSLVTLNLI LAVMAYDRYVATCCPLHYVTAMSPGL CVLLLSLCWGLSVLYGLLLTFLLTRVTFCGPREIHYLFCDMYILLWLACSNTHIIHTAL IATGCFIFLTP LGFMTTSYVRIVRTILQMPSASKKYKTFSTCASHLGVVSLFYGTLAMVYLQPLHTYSMKDSVATVMYAVL TPMMNPFYRLRNKDMHGAPGRVLWRPFQRPK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014566
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014566.1, NP_055381.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	8386
UniProt ID:	<u>P58170</u>
Cytogenetics:	17p13.3
MW:	35.2 kDa