

Product datasheet for **RC224885LIV**

OSTF1 (NM_012383) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	OSTF1
Synonyms:	bA235O14.1; OSF; SH3P2
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_012383
ORF Size:	642 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224885).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_012383.3
RefSeq Size:	901 bp
RefSeq ORF:	645 bp
Locus ID:	26578
UniProt ID:	Q92882
Cytogenetics:	9q21.13



Domains: SH3, ANK

MW: 23.6 kDa

Gene Summary: Osteoclast-stimulating factor-1 is an intracellular protein produced by osteoclasts that indirectly induces osteoclast formation and bone resorption (Reddy et al., 1998 [PubMed 10092216]).[supplied by OMIM, Mar 2008]