

Product datasheet for **RC224889L3V**

NIPAL2 (NM_024759) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NIPAL2
Synonyms:	NPAL2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024759
ORF Size:	1104 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224889).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_024759.1
RefSeq Size:	2255 bp
RefSeq ORF:	1107 bp
Locus ID:	79815
UniProt ID:	Q9H841
Cytogenetics:	8q22.2



Protein Families: Transmembrane

MW: 40.6 kDa