

Product datasheet for **RC224918**

FAM149A (NM_001006655) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FAM149A
Synonyms:	MST119; MSTP119
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC224918 representing NM_001006655
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTATTTGAAGGGAAGTGAACCTCAGACCCAGAGTCTGCTGGCCGAATGCGGGGAGTGACAAGAA
 GATCCCTCCATTTGAGAGTATTAGGAAGACAGCTGATCCTGCCCACTGACAAAGGCGTCCAGCATTTCCA
 GGGCAGCACTCCTGCCTCCGAGTCCACAGACCCCGCTCAGTGCCTGCGGACACAGCAGCAACATCAGA
 GAGTTGTGCATTTCTGGCTCTCAAAATAGTCCCAGCAGCACTCTCAGCCTCTGCCCTGCCAGGCCCTGATG
 ACACAGGGGTTGCTGACCTAACGGCACGTTTCATCCCTGGAAGAAGAGGTTTACCATGTGGATGGAAGAT
 TGAGGAGTATTTTCGCTTTTGACAGAAAAGAGGATGATGATGAATGTCTTGAACAAAAACCACTCAGCCC
 GGTAGGAAATGGCGAACTCGGACTTCTCCTGTTTCCCGCGTGACTGTGTCAAAGATGCCGTGGCAG
 CAGAAGTGTTCGATCAGCTCTGGACAAATATGGTAGAAGTCTTGAAGAGCTGATTAGAAAACTGGGA
 AACTACACTCACAGAAGGGAACAGAGAGAAACATTGAAAGTGGCTGGAACAGATTTCCGCACGTC
 CTCGTTCCACACGCTCAGCCGATGGAGCCAGTGGCCCCCGTCCGGACACGCCGAGGCTCACGGCATCT
 CCCTGGCTTCTCGTCTGAACCCGCCCCAGATTCATCACTTCTCCAGCAGTTTTTATAGTGACATGAATGG
 TGTATGACAATTAAGCAAAACCGCTTCAGCGGAGACCTGCCTATTTTGTGACAGAACGAGAATGAG
 AAGGAGGACAAAGCATCGGGTGGAGGGGAGGTGCTCTCTCCTCCGACCGCACAGACTGGGACGGGCT
 CAGACACTCATGGATTATCACCTTCTGCAAGAAAACACCACTGCCCTGGAGGCTGCCTTCTCTTGTCTT
 AGATTCACAGAGACTAAAACTCCCAACATCTATAGTGACGAAGTCTTCTGGGGAAACAACTGCCGACT
 GGCGTGGACCACATGGCTTCCCACTGGTTCAAACGTCACGGAGCAGGTTCCCCCGCTAGTCACGGAGA
 CCAGGGGGCAGAATACAGCAGTTCCTGGATGCCGCTTGTCTTACAGAGGAAGGCATCTACAAAACCG
 TGTGTTGAGTGCCATGCCTGACGGTACAGAACGATCGCGTCTTCGAGAAAGAACGCCACCTGGAACGG
 TTGTCAAGGCCAGCACAAACCACAGTTCGGCAGTCGGATACGCCTCGAAAAAGTTCATTGACACAAA
 TGGAATTTGCTGCTCACACATGGACAGGTCAAAGTATTTGACAGGTTCAATATGTGCCTAAATCTTT
 TCAGAGGACAATTTGACTTTCAAGAGGAGATTCCAAGTGACATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224918 representing NM_001006655
 Red=Cloning site Green=Tags(s)


MLFEGKVPQTQSLLAECGEWTRRSLHLRVLGRQLILPTDKGVQHFQGSTPASAVHRPPLSACGHSSNIR
 ELCSISGSQIVPAALASALPGPDDTGVADLTARSSLEEEVYHVDGKIEEYFAFDRKEDDDECLEQKPAQP
 GRKWRKLGPPVSPRDCVKDAVAEEVFDHVNTNMVELLEELIRKHWETTLTEGKKQRETLKVAGNRFPHV
 LVPHAHADGASGPPSGHAEAHGISLASRLNPPQIHFFSSSFYSMDMNGVMTIQAKPLQRRPAYFADRTQNE
 KEDKASGGGAGALSSAPHLGRASDTHGLSPSAKKTVPWRLPSLASDSQRLKTPNIYSDEVLRGTKLPT
 GVDHMASPLVQTSRFRPPLVTETRGQNTAVPGCRLVSYRGRHLQNRVLSAMPDGTERTSRLRERTATLER
 LSRPSTHTFRQSDTPRKSSLTQMEFAAHTWTGQSILTGSQYVPKSFQRTTLTFKRRFQVTS

TRTRPLE**QKLI**SEED**LA**AND**ILDYK**DDDD**KV**

Chromatograms: https://cdn.origene.com/chromatograms/mk8018_h03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001006655

ORF Size: 1446 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001006655.3](#)

RefSeq Size: 2422 bp

RefSeq ORF: 1449 bp


Locus ID: 25854

UniProt ID: A5PLN7

Cytogenetics: 4q35.1

MW: 53.3 kDa

Product images:



Circular map for RC224918