

Product datasheet for **RC224924L2V**

Pannexin 3 (PANX3) (NM_052959) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pannexin 3 (PANX3) (NM_052959) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Pannexin 3
Synonyms:	PX3
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_052959
ORF Size:	1176 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224924).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052959.2
RefSeq Size:	1609 bp
RefSeq ORF:	1179 bp
Locus ID:	116337
UniProt ID:	Q96QZ0
Cytogenetics:	11q24.2
Protein Families:	Transmembrane
MW:	44.5 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the innexin family. Innexin family members are known to be the structural components of gap junctions. [provided by RefSeq, Jul 2008]