

Product datasheet for **RC224938A1V**

Human Bim (BCL2L11) (NM_207002) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Bim (BCL2L11) (NM_207002) AAV Particle
Tag: Myc-DDK
Symbol: Bim
Synonyms: BAM; BIM; BOD
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC224938 representing NM_207002
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC
GCC

ATGGCAAAGCAACCTTCTGATGTAAGTTCTGAGTGTGACCGAGAAGGTAGACAATTGCAGCCTGCGGAGA
GGCCTCCCCAGCTCAGACCTGGGGCCCCCTACCTCCCTACAGACAGAGCCACAAGACAGGAGCCCAGCACC
CATGAGTTGTGACAAATCAACACAAACCCCAAGTCCTCCTTGCCAGGCCCTCAACCACTATCTCAGTGCA
ATGGTAGTCATCCTAGAGGATATAGGTGATCTTCACTGTGCTTTGGATTATATTTACTGGCTTAGATT
TGTATGGCCACCACCATAGTCAAGATACAGAACAACCACTCAACCACAAGGATTTCTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224938 representing NM_207002
 Red=Cloning site Green=Tags(s)

MAKQPSDVSSECDREGRQLQPAERPPQLRPGAPTSLQTEPQDRSPAPMSCDKSTQTPSPPCQAFNHLYLSA
MVILEDIGDLSLFCFGFIFTGLDLYGHHSQDTEQLNHKDFS

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_207002
ORF Size: 336 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_207002.1</u>
RefSeq Size:	453 bp
RefSeq ORF:	339 bp
Locus ID:	10018
UniProt ID:	<u>O43521</u>
Cytogenetics:	2q13
MW:	12.2 kDa