

Product datasheet for **RC224948L3V**

ABCB8 (NM_007188) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ABCB8 (NM_007188) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABCB8
Synonyms:	EST328128; M-ABC1; MABC1; MITOSUR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007188
ORF Size:	2154 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224948).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007188.2 , NP_009119.1
RefSeq Size:	2439 bp
RefSeq ORF:	2157 bp
Locus ID:	11194
UniProt ID:	Q9NUT2
Cytogenetics:	7q36.1
Domains:	ABC_membrane, ABC_tran, AAA
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: ABC transporters

MW: 77.9 kDa

Gene Summary: This nuclear gene encodes a multi-pass membrane protein that is targeted to the mitochondrial inner membrane. The encoded protein is an ATP-dependent transporter that may mediate the passage of organic and inorganic molecules out of the mitochondria. Loss of function of the related gene in mouse results in a disruption of iron homeostasis between the mitochondria and cytosol. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]