

## Product datasheet for **RC224960**

### Elastin (ELN) (NM\_001081755) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Elastin (ELN) (NM_001081755) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ELN
Synonyms:	ADCL1; SVAS; WBS; WS
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RC224960 representing NM\_001081755  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGCGGGTCTGACGGCGCGGCCCGCGGCCGGAGTCTCCTGCTCCTGCTGTCCATCTCCACCCT  
CTCGGCCTGGAGGGTCCCTGGGGCCATTCTGGTGGAGTTCCTGGAGGAGTCTTTTATCCAGGGCTGG  
TCTCGGAGCCCTTGAGGAGGAGCGCTGGGGCCTGGAGGCAAACCTTTAAGCCAGTTCGCGAGGGCTT  
GCGGGTGTGGCCTTGGGGCAGGGCTCGGGCCTTCCCGCAGTTACCTTCCGGGGCTCTGGTGCCTG  
GTGGAGTGGCTGACGCTGCTGCAGCCTATAAAGCTGCTAAGGCTGGCGCTGGGCTTGGTGGTGTCCAGG  
AGTTGGTGGCTTAGGAGTGTCTGACGTTGCGGTGGTTCCTCAGCCTGGAGCCGGAGTGAAGCCTGGGAAA  
GTGCCGGTGTGGGGCTGCCAGGTGTATACCCAGGTGGCGTGTCCAGGAGCTCGGTTCCCGGTGTGG  
GGGTGCTCCCTGGAGTTCCTACTGGAGCAGGAGTTAAGCCAAAGCTCCAGGTGTAGGTGGAGCTTTTGC  
TGGAAATCCAGGAGTTGGACCTTTGGGGACCGCAACCTGGAGTCCACTGGGGTATCCCATCAAGGCC  
CCCAAGCTGCCTGGTGGCTATGGACTGCCTACACCACAGGGAACTGCCCTATGGCTATGGGCCCGGAG  
GAGTGGCTGGTGCAGCGGCAAGGCTGGTTACCCAACAGGGACAGGGGTTGGCCCCAGGCAGCAGCAGC  
AGCGGCAGCTAAAGCAGCAGCAAAGTTCGGTGTGGAGCAGCCGGAGTCTCCCTGGTGTGGAGGGCT  
GGTGTTCCTGGCGTGCCTGGGGCAATTCCTGGAATTGGAGGATCGCAGGCGTTGGGACTCCAGCTGCAG  
CTGCAGCTGCAGCAGCAGCCGCTAAGGCAGCAAGTATGGAGTGTCTGCAGGCTTAGTGCCTGGTGGCC  
AGGCTTTGGCCGGGAGTAGTTGGTGTCCAGGAGCTGGCGTTCAGGTGTTGGTGTCCAGGAGCTGGG  
ATTCAGTGTCCAGGTGTGGGATCCAGGTGTGGGTTCCAGGGTGTGTACCAGAAAGCAGCTG  
CTAAGGCAGCTGCAAAGGCAGCCAAATACGGGGCCAGGCCGGAGTCGGAGTTGGAGGCACTTCTACTTA  
CGGGGTTGGAGCTGGGGCTTTCCCGGCTTTGGTGTGGAGTCCGAGGATCCCTGGAGTGCAGGTGTC  
CCTGGTGTGGAGGTGTTCCCGAGTCCGAGGTGTCCCGGAGTTGGCATTTCGCCGAAGCTCAGGCAG  
CAGCTGCCGCAAGGCTGCCAAGTACGGGTTAGTTCCTGGTGTGGCGTGGCTCCTGGAGTTGGCGTGGC  
TCCTGGTGTGGTGTGGCTCCTGGAGTTGGCTTGGCTCCTGGAGTTGGCGTGGCTCCTGGAGTTGGTGTG  
GCTCCTGGCGTTGGCGTGGCTCCCGCATTGGCCCTGGTGGAGTGCAGCTGCAGCAAAATCCGCTGCCA  
AGGTGGCTGCCAAAGCCAGCTCCGAGCTGCAGTGGGCTTGGTGTGGCATCCCTGGACTTGGAGTTGG  
TGTCCGGCTCCCTGGACTTGGAGTTGGTGTGGTTCCTGGACTTGGAGTTGGTGTGGTGTTCCTGGC  
TTCGGGGCAGTACCTGGAGCCCTGGCTGCCGCTAAAGCAGCCAAATATGGAGCAGCAGTGCCTGGGTCC  
TTGGAGGGCTCGGGCTCTCGGTGGAGTAGGCATCCAGGCGGTGTGGTGGGAGCCGGACCCGCCCGCC  
CGCTGCCGAGCCAAAGCTGCTGCCAAAGCCGCCAGTTTGGCCTAGTGGGAGCCGCTGGGCTCGGAGGA  
CTCGGAGTCGGAGGGCTTGGAGTTCAGGTGTTGGGGCCTTGGAGGTATACCTCCAGCTGCAGCCGCTA  
AAGCAGCTAAATACGGTGTCTGGCCTTGGAGGTGTCTAGGGGTGCCGGCAGTTCCCACTTGGAGG  
AGTGGCAGCAAGACCTGGCTTCGGATTGTCTCCATTTCCAGGTGGGGCTGCCTGGGAAAGCTTGT  
GGCCGGAAGAGAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224960 representing NM\_001081755  
 Red=Cloning site Green=Tags(s)


MAGLTAAAPRPGVLLLLLSILHPSRPGGVPGAIPGGVPGGVFYPGAGLGALGGGALGPGGKPLKPVPGGL  
 AGAGLGAGLGAFPAVTFPGALVPGGVADAAAAYKAAKAGAGLGGVPGVGGGLVVSAGAVVPPQGAGVKPGK  
 VPGVGLPGVYPPGGVLPGARFPVGVVLPVPTGAGVKPKAPGVGGAFAGIPGVGPFGGPQPGVPLGYPIKA  
 PKLPGGYGLPYTTGKLPYGYGPGGVAGAAGKAGYPTGTGVGPQAAAAAAKAAAKFGAGAAGVLPVGGGA  
 GVPVPGAIPGIIGGIAGVGTAAAAAAYKAAKAGAGLGGVPGVGGGLVVSAGAVVPPQGAGVKPGK  
 IPVVPAGIPGAAVPGVVSPEAAAAKAAKAAKYGARPGVGGIPTYGVGAGGFPFGVGVGGIPGVAGV  
 PGVGGVPGVGGVPGVVISPEAQAAAAKAAKYGLVPGVGVAPGVGVAPGVGVAPGVGLAPGVGVAPGVGV  
 APGVGVAPGIGPGGVAAAAKSAKVAQAQLRAAAGLGGVGGVGGVGGVGGVGGVGGVGGVGGVGGVGGV  
 FGAVPGALAAKAAKYGAAVPGVGLGGLGALGGVGGVGGVGGVGGVGGVGGVGGVGGVGGVGGVGGVGGV  
 LGVGGV  
 GRKRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



\* The last codon before the Stop codon of the ORF

**Plasmid Map:**


**ACCN:** NM\_001081755

**ORF Size:** 2115 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001081755.3](#)

**RefSeq Size:** 3423 bp

**RefSeq ORF:** 2118 bp

<b>Locus ID:</b>	2006
<b>UniProt ID:</b>	<a href="#">P15502</a>
<b>Cytogenetics:</b>	7q11.23
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>MW:</b>	60.98 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein that is one of the two components of elastic fibers. Elastic fibers comprise part of the extracellular matrix and confer elasticity to organs and tissues including the heart, skin, lungs, ligaments, and blood vessels. The encoded protein is rich in hydrophobic amino acids such as glycine and proline, which form mobile hydrophobic regions bounded by crosslinks between lysine residues. Degradation products of the encoded protein, known as elastin-derived peptides or elastokines, bind the elastin receptor complex and other receptors and stimulate migration and proliferation of monocytes and skin fibroblasts. Elastokines can also contribute to cancer progression. Deletions and mutations in this gene are associated with supravalvular aortic stenosis (SVAS) and autosomal dominant cutis laxa. [provided by RefSeq, Aug 2017]</p>