

## Product datasheet for **RC224976L2V**

### USP17L2 (NM\_201402) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	USP17L2 (NM_201402) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP17L2
Synonyms:	DUB-3; DUB3; USP17
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_201402
ORF Size:	1590 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224976).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_201402.1</a>
RefSeq Size:	1593 bp
RefSeq ORF:	1593 bp
Locus ID:	377630
UniProt ID:	<a href="#">Q6R6M4</a>
Cytogenetics:	8p23.1
MW:	59.4 kDa



[View online »](#)

**Gene Summary:**

DUB3 is a member of the ubiquitin processing protease (UBP) subfamily of deubiquitinating enzymes. See USP1 (MIM 603478) for background information.[supplied by OMIM, Mar 2008]