

## Product datasheet for **RC225001L4V**

### serum amyloid A2 (SAA2) (NM\_001127380) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	serum amyloid A2 (SAA2) (NM_001127380) Human Tagged ORF Clone Lentiviral Particle
Symbol:	serum amyloid A2
Synonyms:	SAA; SAA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001127380
ORF Size:	249 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225001).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001127380.2</a>
RefSeq Size:	2053 bp
RefSeq ORF:	252 bp
Locus ID:	6289
UniProt ID:	<a href="#">P02735</a>
Cytogenetics:	11p15.1
Protein Families:	Druggable Genome
MW:	9.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the serum amyloid A family of apolipoproteins. The encoded preproprotein is proteolytically processed to generate the mature protein. This protein is a major acute phase protein that is highly expressed in response to inflammation and tissue injury. This protein also plays an important role in HDL metabolism and cholesterol homeostasis. High levels of this protein are associated with chronic inflammatory diseases including atherosclerosis, rheumatoid arthritis, Alzheimer's disease and Crohn's disease. This protein may also be a potential biomarker for certain tumors. Finally, antimicrobial activity against *S. aureus* and *E. coli* resides in the N-terminal portion of the mature protein. [provided by RefSeq, Jul 2020]