

Product datasheet for **RC225009L3V**

ZNF673 (KRBOX4) (NM_001129900) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | ZNF673 (KRBOX4) (NM_001129900) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ZNF673 |
| Synonyms: | ZNF673 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001129900 |
| ORF Size: | 279 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC225009). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001129900.1 , NP_001123372.1 |
| RefSeq ORF: | 282 bp |
| Locus ID: | 55634 |
| UniProt ID: | Q5JUW0 |
| Cytogenetics: | Xp11.3 |
| Protein Families: | Transcription Factors |
| MW: | 10.3 kDa |



[View online »](#)

Gene Summary:

This encodes a zinc finger protein with an N-terminal KRAB (Krüppel-associated) domain found in transcriptional repressors. This gene is located in a region of the X chromosome thought to be involved in nonsyndromic X-linked cognitive disability. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]