

Product datasheet for **RC225018A1V**

Human LINC02877 (NM_001101337) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LINC02877 (NM_001101337) AAV Particle
Tag:	Myc-DDK
Symbol:	LINC02877
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC225018 representing NM_001101337 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC ATGCATGTATGTGATCTGCAGAAGCTGGTCAGGATCCAACCTTGCATTTACTACATTTCCATGGATGTTCA GTTGCCATCTTCTGCCTACTCCTGAGCTTTCATCCAAGAGAAATCAGTGTGTTGCTTTACAAGACTTCCGG TTGTCTAACCCAAATGCCAATTTTATATGGGCACCCTGCCACTTCTAAAAGACTACATTTGCAGGCA ATTTTGAACCCAGGAAAGAAGATTCAAGGAGGAACTGAAATACAGAGAGGCTCCTTTGCAAACCCAGTATC AAACAGATGCCAGTCACTTG ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC225018 representing NM_001101337 Red=Cloning site Green=Tags(s) MHVCDLQKLVRIQLAFTTFPMFSCHELLPTPELSSKRNQCLLYKTSGLTQMPILYGHPTLLKDYILQA ILQPGKKIQGGTEIQRGSFANQYQTDASHL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001101337
ORF Size:	300 bp
Buffer:	PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_001101337.1](#), [NP_001094807.1](#)

RefSeq ORF: 303 bp

Locus ID: 152118

Cytogenetics: 3q25.2

MW: 11.2 kDa