

Product datasheet for **RC225060A1V**

Human YPEL5 (NM_001127400) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human YPEL5 (NM_001127400) AAV Particle
Tag: Myc-DDK
Symbol: YPEL5
Synonyms: CGI-127
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225060 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGCAGAATTTTCCTTGATCATATCGGTGGTACCCGTCGTTTTCTTGTGCAAACGTGATACGATCC
TGACCAACCGCTCAGAACTCATCTCCACTCGTTTCACAGGCGCCACTGGCAGAGCATTTCTTTTAA
GGTAGTTAACCTGCAGTACAGTGAAGTTCAAGATCGGGTCATGCTCACTGGCCGCCACATGGTTCGAGAT
GTGAGCTGCAAACTGCAATAGCAAACCTGGGATGGATCTATGAGTTGCCACTGAAGACAGCCAGCGAT
ATAAGGAAGCCCGGTGATCCTGGAACGTGCTCTAGTTCGAGAGAGTGAGGGCTTTGAGGAGCATGTACC
ATCTGATAACTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225060 protein sequence
Red=Cloning site Green=Tags(s)

MGRIFLDHIGGTRLFSCANCDTILTNRSELISTRFTGATGRAFLFNKVVNLQYSEVQDRVMLTGRHMVRD
VSKKNCNSKLGWIYEFATEDSQRYKEGRVILERALVRESEGFEHVPDMS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001127400



[View online »](#)

ORF Size:	363 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001127400.1</u>
RefSeq Size:	2578 bp
RefSeq ORF:	366 bp
Locus ID:	51646
UniProt ID:	<u>P62699</u>
Cytogenetics:	2p23.1
MW:	13.8 kDa