

Product datasheet for RC225061A1V

Human YPEL5 (NM_001127401) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human YPEL5 (NM_001127401) AAV Particle
Tag: Myc-DDK
Symbol: YPEL5
Synonyms: CGI-127
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225061 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCAGAATTTTCCTTGATCATATCGGTGGTACCCGCTGTTTTCTTGTGCAAACGTGATACGATCC
 TGACCAACCGCTCAGAACTCATCTCCACTCGTTTCACAGGCGCCACTGGCAGAGCATTTCTTTTAA
 GGTAGTTAACCTGCAGTACAGCGAAGTTCAAGATCGGGTCATGCTCACTGGCCGCCACATGGTTCGAGAT
 GTGAGCTGCAAAACTGCAATAGCAAACCTGGGATGGATCTATGAGTTGCCACTGAAGACAGCCAGCGAT
 ATAAGGAAGCCCGGTGATCCTGGAACGTGCTCTAGTTCGAGAGAGTGAGGGCTTTGAGGAGCATGTACC
 ATCTGATAACTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225061 protein sequence
 Red=Cloning site Green=Tags(s)

MGRIFLDHIGGTRLFSCANCDTILTNRSELISTRFTGATGRAFLFNKVVNLQYSEVQDRVMLTGRHMVRD
 VSKKNCNSKLGWIYEFATEDSQRYKEGRVILERALVRESEGFEHVPDSDNS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001127401



[View online »](#)

ORF Size:	363 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001127401.1</u> , <u>NP_001120873.1</u>
RefSeq Size:	2639 bp
RefSeq ORF:	366 bp
Locus ID:	51646
UniProt ID:	<u>P62699</u>
Cytogenetics:	2p23.1
MW:	13.8 kDa