

Product datasheet for **RC225070L3V**

PSMG4 (NM_001128591) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PSMG4 (NM_001128591) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PSMG4
Synonyms:	bA506K6.2; C6orf86; PAC4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001128591
ORF Size:	369 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225070).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001128591.1
RefSeq ORF:	372 bp
Locus ID:	389362
UniProt ID:	Q5JS54
Cytogenetics:	6p25.2
MW:	13.6 kDa
Gene Summary:	Chaperone protein which promotes assembly of the 20S proteasome.[UniProtKB/Swiss-Prot Function]


[View online »](#)