

Product datasheet for **RC225088L1V**

KCNE1 (NM_001127670) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	KCNE1
Synonyms:	ISK; JLNS; JLNS2; LQT2/5; LQT5; Mink
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001127670
ORF Size:	387 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225088).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001127670.2
RefSeq Size:	3424 bp
RefSeq ORF:	390 bp
Locus ID:	3753
UniProt ID:	P15382
Cytogenetics:	21q22.12



Protein Families: Druggable Genome, Ion Channels: Other, Transmembrane

MW: 15.1 kDa

Gene Summary: The product of this gene belongs to the potassium channel KCNE family. Potassium ion channels are essential to many cellular functions and show a high degree of diversity, varying in their electrophysiologic and pharmacologic properties. This gene encodes a transmembrane protein known to associate with the product of the KVLQT1 gene to form the delayed rectifier potassium channel. Mutation in this gene are associated with both Jervell and Lange-Nielsen and Romano-Ward forms of long-QT syndrome. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]