

Product datasheet for **RC225108A1V**

Human PIRT (NM_001101387) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PIRT (NM_001101387) AAV Particle
Tag: Myc-DDK
Symbol: PIRT
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225108 representing NM_001101387
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACGATGGAGACTCTCCCAAGGTTCTAGAGGTCGATGAGAAGTCTCCAGAAGCCAAGGACCTGCTGC
 CCAGCCAGACCGCCAGCTCCCTGTGCATCAGCTCCAGGAGCGAGTCTGTCTGGACCACCACCCAGGAG
 TAACTGGGAAATCTACCGCAAGCCCATCGTTATCATGTCAGTGGGCGGTGCCATCCTGCTTTTCGGCGTG
 GTCATCACCTGCTTGGCCTACACCTTGAAGCTGAGTGACAAGAGTCTCTCCATCCTCAAATGGTAGGGC
 CTGGCTTCTGTCCCTGGGACTCATGATGCTGGTGTGCGGGCTGGTGTGGGTGCCATCATCAAAAAGAA
 ACAGAAGCACAGACAGAAGTCGAATTTCTACGCAGCCTCAAGTCTTCTCTGACTCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225108 representing NM_001101387
 Red=Cloning site Green=Tags(s)

MTMETLPKYLEVDEKSPEAKDLLPSQTASSLCISSRSESVWTTTPRSNWEIYRKPVIMSVGGAILLFGV
 VITCLAYTLKLSDKSLSILKMGVPGFSLGLMMLVCGLVWVPIIKKKQKHRQKSNFLRSLKSFLLTR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001101387
ORF Size: 411 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001101387.1</u>
RefSeq ORF:	414 bp
Locus ID:	644139
UniProt ID:	<u>POC851</u>
Cytogenetics:	17p12
MW:	15.2 kDa