

Product datasheet for **RC225122L3V**

FXYD2 (NM_001127489) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | FXYD2 (NM_001127489) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | FXYD2 |
| Synonyms: | ATP1G1; HOMG2; MGC12372 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001127489 |
| ORF Size: | 435 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC225122). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001127489.1 , NP_001120961.1 |
| RefSeq ORF: | 437 bp |
| Locus ID: | 486 |
| Cytogenetics: | 11q23.3 |
| Protein Families: | Druggable Genome, Ion Channels: Other, Transmembrane |
| MW: | 14.9 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the FXVD family of transmembrane proteins. This particular protein encodes the sodium/potassium-transporting ATPase subunit gamma. Mutations in this gene have been associated with Renal Hypomagnesemia-2. Alternatively spliced transcript variants have been described. Read-through transcripts have been observed between this locus and the upstream FXVD domain-containing ion transport regulator 6 (FXVD6, GeneID 53826) locus.[provided by RefSeq, Feb 2011]