

Product datasheet for **RC225126L4V**

FAM72A (NM_001123168) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: FAM72A

Synonyms: LMP1P; p17; Ugene

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_001123168

ORF Size: 447 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225126).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001123168.1](#), [NP_001116640.1](#)

RefSeq ORF: 450 bp

Locus ID: 729533

UniProt ID: [Q5TYM5](#)

Cytogenetics: 1q32.1

MW: 16.4 kDa



Gene Summary:

May play a role in the regulation of cellular reactive oxygen species metabolism. May participate in cell growth regulation.[UniProtKB/Swiss-Prot Function]