

Product datasheet for RC225130A1V

Human LINC02902 (NM_001123065) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LINC02902 (NM_001123065) AAV Particle
Tag: Myc-DDK
Symbol: LINC02902
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225130 representing NM_001123065
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAGGATGGCACCAACAGAGTCAACAGAAGGAAGACGGCTCTGGCCGGGCCCCAGGGAGGGAGGCAGCG
 GTAAGGAAACAACCTTCAGAGAAGTTAAGCAACTTGCCAGGCCACACAGCTATTCACCAAAGAGAGCTGA
 TGCTGAGTCTTTCAGAGGAGTGCCTGCAGCATTTAAAAATGCAGAGAAGTGTTCAGAGCCTGCTGGGGA
 AGCAGGGAGCTGCTATTTCTGTTCAAGCAATCAGTGAGGCTGGACCTGCCAGAATTCATGTGGAATCA
 CCCTAGAGAAGGCTGGTGGCTTGAAGACACTGGATCTCACTGGCTCAGCTGGGCACGGTGCAAGGTGCT
 ATACATAAATGGTTTCACTGACCCTGGAAGGATGCTCAGGCCTGGATACTCATTGTGAGCTGCAAAAAA
 GGAAAGGGGACCCTGAGAGGGAAGGCAGGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225130 representing NM_001123065
 Red=Cloning site Green=Tags(s)

MRMAPTESTEGRRLWPGPREGGSGKETTSEKLSNLPRPHSYSPKRADAESFRGVPAAFKKCREVFRACWG
 SRELLFLFKAISEAGPAQNSCGITLEKAGGLEDTGSHLWSWARCKVLYINGFTDPWKDAQAWILIVSCKK
 GKGTPEREGRN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001123065


[View online »](#)

ORF Size:	453 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001123065.1</u> , <u>NP_001116537.1</u>
RefSeq ORF:	456 bp
Locus ID:	401335
UniProt ID:	<u>Q6ZTY9</u>
Cytogenetics:	7p12.3
MW:	16.6 kDa