

Product datasheet for **RC225132L3V**

TMEM166 (EVA1A) (NM_001135032) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: TMEM166

Synonyms: FAM176A; TMEM166

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_001135032

ORF Size: 456 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC225132).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001135032.1](#)

RefSeq Size: 1872 bp

RefSeq ORF: 459 bp

Locus ID: 84141

UniProt ID: [Q9H8M9](#)

Cytogenetics: 2p12



Protein Families: Transmembrane

MW: 17.5 kDa

Gene Summary: Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.
[UniProtKB/Swiss-Prot Function]