

Product datasheet for **RC225138**

RABL5 (IFT22) (NM_001130820) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RABL5 (IFT22) (NM_001130820) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RABL5
Synonyms: FAP9; RABL5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC225138 representing NM_001130820
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCTGAAAGCCAAGATCCTCTTCGTGGGGCCTTGCGAGAGTGGAAAACTGTTTTGGCCAACTTTCTGA
 CAGAATCTTCTGACATCACTGAATACAGCCCAACCAAGGAGTGAGGTTTGAGTCCTGCTGGCCGGCCCT
 GATGAAGGATGCTCATGGAGTGGTGATCGTCTTCAATGCTGACATCCCAAGCCACCGGAAGGAAATGGAG
 ATGTGGTATTCTGCTTTGTCCAACAGCCGTCCTTACAGGACACACAGTGTATGCTAATTGCACACCACA
 AACCAGGCTCTGGAGATGATAAAGGAAGCCTGTCTTTGTCGCCACCCTTGAACAAGCTGAAGCTGGTGCA
 CTCAAACCTGGAAGATGACCCTGAGGAGATCCGGATGGAATTCATAAAGTATTTAAAAAGCATAATCAAC
 TCCATGTCTGAGAGCAGAGACAGGGAGGAGATGTCAATTATGACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC225138 representing NM_001130820
Red=Cloning site Green=Tags(s)

MLKAKILFVGPCESGKTVLANFLTESSDITEYSPTQGVRFESCWPALMKDAHGVVIVFNADIPSHRKEME
 MWYSCFVQQPSLQDTQCMLIAHHKPGSGDDKGSLSLSPPLNKLKVHNSLEDDPEEIRMEFIKYLKSIIN
 SMSES RDREMSIMT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_001130820

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001130820.3](#)


RefSeq ORF: 468 bp

Locus ID: 64792

UniProt ID: [Q9H7X7](#)

Cytogenetics:	7q22.1
Protein Families:	Druggable Genome
MW:	17.4 kDa
Gene Summary:	Small GTPase-like component of the intraflagellar transport (IFT) complex B.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC225138