

Product datasheet for **RC225167A1V**

Human CNBP (NM_001127196) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CNBP (NM_001127196) AAV Particle
Tag:	Myc-DDK
Symbol:	CNBP
Synonyms:	CNBP1; DM2; PROMM; RNF163; ZCCHC22; ZNF9
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC225167 representing NM_001127196 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGCAGCAATGAGTGCTTCAAGTGTGGACGATCTGGCCACTGGGCCCGGAATGTCCTACTGGTGGAG
GCCGTGGTCGTGGAATGAGAAGCCGTGGCAGAGGTTCCAGTTTGTTCCTCGTCTCTCCAGACATTTG
TTATCGCTGTGGTGAGTCTGGTCATCTTGCCAAGGATTGTGATCTTCAGGAGGATGCCTGCTATAACTGC
GGTAGAGGTGGCCACATTGCCAAGGACTGCAAGGAGCCCAAGAGAGAGCGAGAGCAATGCTGCTACAAC
GTGGCAAACAGGCCATCTGGCTCGTGACTGCGACCATGCAGATGAGCAGAAATGCTATTCTTGTGGAGA
ATTCGGACACATTCAAAAAGACTGCACCAAAGTGAAGTGTATAGGTGTGGTGAAGTGGTCAATGTAGCC
ATCAACTGCAGCAAGACAAGTGAAGTCAACTGTTACCGCTGTGGCGAGTCAGGGCACCTTGCACGGGAAT
GCACAATTGAGGCTACAGCC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC225167 representing NM_001127196 Red=Cloning site Green=Tags(s)
-------------------	---

MSSNECFKCGRSGHWARECPTGGGRGRGMRSRGRGFQFVSSSLPDICYRCGESGHLAKDCDLQEDACYNC
GRGGHIAKDCKEPREREQCCYNCCKPGHLARCDHADEQKCYSCGEFGHIQKDKTKVKCYRCGETGHVA
INCSKTSEVNCYRCGESGHLARECTIEATA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species:	Human
Serotype:	AAV-2



[View online »](#)

ACCN:	NM_001127196
ORF Size:	510 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001127196.1</u>
RefSeq ORF:	513 bp
Locus ID:	7555
UniProt ID:	<u>P62633</u>
Cytogenetics:	3q21.3
MW:	18.6 kDa