

Product datasheet for **RC225181L3V**

JDP2 (NM_001135049) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	JDP2 (NM_001135049) Human Tagged ORF Clone Lentiviral Particle
Symbol:	JDP2
Synonyms:	JUNDM2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001135049
ORF Size:	522 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225181).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001135049.1 , NP_001128521.1
RefSeq Size:	3960 bp
RefSeq ORF:	525 bp
Locus ID:	122953
UniProt ID:	Q8WYK2
Cytogenetics:	14q24.3
Protein Families:	Transcription Factors
MW:	19.8 kDa



[View online »](#)

Gene Summary:

Component of the AP-1 transcription factor that represses transactivation mediated by the Jun family of proteins. Involved in a variety of transcriptional responses associated with AP-1 such as UV-induced apoptosis, cell differentiation, tumorigenesis and antitumorigenesis. Can also function as a repressor by recruiting histone deacetylase 3/HDAC3 to the promoter region of JUN. May control transcription via direct regulation of the modification of histones and the assembly of chromatin.[UniProtKB/Swiss-Prot Function]