

Product datasheet for RC225182A1V

Human CLEC2A (NM_001130711) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CLEC2A (NM_001130711) AAV Particle
Tag: Myc-DDK
Symbol: CLEC2A
Synonyms: INPE5792; KACL; PILAR; UNQ5792
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225182 representing NM_001130711
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGATTAATCCAGAGCTGCGGGATGGCAGAGCTGATGGCTTCATACATCGGATAGTTCCCAAGTTGATAC
 AAACTGGAAGATTGGCCTTATGTGCTTCCTGAGTATTATTACTACAGTTTGCAATTATTATGATAGC
 CACATGGTCCAAGCATGCTAAACCTGTGGCATGTTCAAGGGACTGGCTTGGAGTGAGAGATAAGTGTTTC
 TATTTTCTGATGATACCAAGAAATTGGACAGCCAGTAAATATTTGTAGTTTGAGAAAGCAGAACTTG
 CTCAGATTGATACACAAGAAGACATGGAATTTTGAAGAGGTACGCAGGAAGTATGCACTGGATTGG
 ACTAAGCAGGAAACAAGGAGATTCTTGAAATGGACAAATGGCACCACATTCAATGGTTGGTTTGAAATT
 ATAGGGAACGGATCCTTTGCTTTCTTGAGTGCTGATGGAGTCCATAGTCCAGAGGATTTATTGATATCA
 AGTGGATTTGCAGCAACCTAAATATTTTAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225182 representing NM_001130711
 Red=Cloning site Green=Tags(s)

MINPELRDGRADGFIHRIVPKLIQNWIGLMCFLSIIITTVCIIMIATWSKHAKPVACSGDWLGVRDKCF
 YFSDDRNWTASKIFCSLQKAELAQIDTQEDMEFLKRYAGTDMHWIGLSRKQGDGSKWTNGTTFNWFEI
 IGNGSFAFLSADGVHSSRGFIDIKWICKPKYFL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_001130711
ORF Size:	522 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001130711.1</u>
RefSeq ORF:	525 bp
Locus ID:	387836
UniProt ID:	<u>Q6UVW9</u>
Cytogenetics:	12p13.31
MW:	19.8 kDa