

Product datasheet for **RC225182L4V**

CLEC2A (NM_001130711) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CLEC2A (NM_001130711) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLEC2A
Synonyms:	INPE5792; KACL; PILAR; UNQ5792
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001130711
ORF Size:	522 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC225182).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130711.1
RefSeq ORF:	525 bp
Locus ID:	387836
UniProt ID:	Q6UWV9
Cytogenetics:	12p13.31
MW:	19.8 kDa
Gene Summary:	CLEC2A belongs to the CLEC2 family of activation-induced, natural killer gene complex-encoded C-type lectin-like receptors (Spreu et al., 2007 [PubMed 18046548]).[supplied by OMIM, May 2008]


[View online »](#)