

Product datasheet for **RC225204A1V**

Human PHR1 (PLEKHB1) (NM_001130035) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PHR1 (PLEKHB1) (NM_001130035) AAV Particle
Tag: Myc-DDK
Symbol: PHR1
Synonyms: KPL1; PHR1; PHRET1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC225204 representing NM_001130035
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCCTGGTGAGGGGCGGCTGGCTGTGGAGACAGAGCTCCATCCTCCGCCGCTGGAAGCGGAAGTGGT
TTGCCCTGTGGCTGGACGGGACCCTGGGATACTACCACGATGAGACAGCGCAGGACGAGGAGGACCGTGT
GCTCATCCACTTCAATGTCCGTGACATAAAGATCGGCCAGAGTGCCATGATGTGCAGCCCCAGAGGGC
CGGAGCCGAGATGGCTGCTGACTGTGAACCTACGGGAAGGCGGCCGCTGCACCTCTGTGCGGAGACCA
AGGATGATGCCCTAGCATGGAAGACAGCACTGCTGGAGGCAAACTCCACCCCGGTGCGCGTCTACAGCCC
GTACCAAGACTACTACGAGGTGGTGCCCCCAATGCACACGAGGCCACGTATGTCCGCAGCTACTACGGA
CCGCCCTACGACGGCCCTGGCGTGACGCACGTGATAGTGCGGGAGGATCCCTGCTACAGCGCCGGCGCCC
CTCTGGCCATGGGCATGCTTGCGGGAGCCGCACTGGGCGGCGCTGGGCTCGCTCATGTGGTCGCCCTG
CTGGTTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC225204 representing NM_001130035
Red=Cloning site Green=Tags(s)

MALVRGGWLWRQSSILRRWKRNFALWLDGTLGYHDETAQDEEDRVLIHFNVRDIKIGPECHDVQPPEG
RSRDGLLTVNLREGGRLHLCAETKDDALAWKTALLEANSTPVRVYSPYQDYEVVPPNAHEATYVRSYYG
PPYAGPGVTHIVREDPCYSAGAPLAMGMLAGAATGAALGSLMWSPCWF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001130035
ORF Size:	567 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001130035.1</u>
RefSeq ORF:	570 bp
Locus ID:	58473
UniProt ID:	<u>Q9UF11</u>
Cytogenetics:	11q13.4
MW:	21.1 kDa